


PicoLITH, Agilent Lithography and Nanomanipulation Package



Data Sheet



Features and Benefits

- CAD-style drawing program offers user-friendly interface for easy operation
- Ability to design and save shapes for reuse helps conserve time and effort
- Ability to control parameters such as force or probe bias increases utility
- Ability to work “on-the-fly” or through scripts provides enhanced flexibility
- Move-tip tool provides real-time lithography
- C/C++, Labview, MATLAB, Python scripting enables custom, programmable lithography
- Seamless integration with Agilent’s AFM systems

Additional Capabilities

- Micro-machining via force variations
- Nano-patterning via probe bias and EC etching
- Nano-writing via surface potential controlled deposition and modification

Overview

Agilent’s lithography and nanomanipulation software is a CAD-style package that provides a nanolithography tool for scanning probe microscopy (SPM) researchers. It allows users to sketch various shapes on a canvas (including lines, poly-lines, circles, and boxes) that can then be mapped to a real sample surface. Individual shapes can be assigned a unique set of parameters. Each parameter set is associated



Agilent Technologies

with a color for easy visual reference. Accessible parameters include force, probe bias, sample bias or potential, and probe speed. As shown in Figure 1, a user-friendly interface enables the user to edit parameters and segment properties, such as start and end points and the order of execution. The software offers users the ability to cut and paste predesigned shapes to assist in the building of more complex patterns. Once designed, shapes can be reused in future experiments and are easily resized.

Move-Tip Tool Applications

The move-tip tool provides a complete set of functions for nanomanipulation. With numerous parameter settings, such as force set-point or tip bias, the user can inscribe a line, inscribe a pattern, or manipulate particles in real time. (See Figures 2 and 3.)




Figure 1. PicoLITH design canvas (top) over a grey-scale AFM image, primitive property window (left), parameter settings window (right).




Figure 2. Making a smiling face with gold particles on a mica surface. Scan size: 2.5 µm x 2.5 µm.




Figure 3. Sequence of drawing triangles on CD blank surface. Green lines indicate the path of the tip movement in contact mode. Images are taken in AAC mode. Scan size: 5.5 µm x 5.5 µm.

Applications

Agilent AFM Systems come with a 32–32 bit Windows XP-based software package that incorporates scripting capabilities. Scripting allows the user to program most of the functions in the microscope software and define special sequences of desired action. By combining Agilent’s

closed-loop scanners and stages with scripting capability, the user is provided with a unique scripting tool for lithography applications. With predefined parameter settings (including STM/AFM tip bias, force set-point, and more), scripting makes drawing a line or a circle on the sample surface a straightforward procedure (Figure 4).




Figure 4a. AAC mode topography image of PZT film. Scan size: 14 μm x 14 μm .




Figure 4b. Using Agilent’s lithography and manipulation software to write with tip bias (left). Drawing and parameter setting (right).




Figure 4c. EFM image after writing. Scan size: 14 μm x 14 μm .




Figure 5. Drawing circles on PMMA film. Scan size: 2.5 μm x 2.5 μm .




Figure 6. Drawing of Langzhou University logo.




Figure 7. Inscribing a circle and lines on a polymer film by increasing the tip-sample interaction force. Scan size: 5 μm x 5 μm .

AFM Instrumentation from Agilent Technologies

Agilent Technologies offers high-precision, modular AFM solutions for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Agilent's leading-edge R&D laboratories are dedicated to the timely introduction and optimization of innovative and easy-to-use AFM technologies.

www.agilent.com/find/afm

Americas

Canada	(877) 894 4414
Latin America	305 269 7500
United States	(800) 829 4444

Asia Pacific

Australia	1 800 629 485
China	800 810 0189
Hong Kong	800 938 693
India	1 800 112 929
Japan	0120 (421) 345
Korea	080 769 0800
Malaysia	1 800 888 848
Singapore	1 800 375 8100
Taiwan	0800 047 866
Thailand	1 800 226 008

Europe & Middle East

Austria	43 (0) 1 360 277 1571
Belgium	32 (0) 2 404 93 40
Denmark	45 70 13 15 15
Finland	358 (0) 10 855 2100
France	0825 010 700*
	*0.125 €/minute
Germany	49 (0) 7031 464 6333
Ireland	1890 924 204
Israel	972-3-9288-504/544
Italy	39 02 92 60 8484
Netherlands	31 (0) 20 547 2111
Spain	34 (91) 631 3300
Sweden	0200-88 22 55
Switzerland	0800 80 53 53
United Kingdom	44 (0) 118 9276201

Other European Countries:

www.agilent.com/find/contactus

Product specifications and descriptions in this document subject to change without notice.

© Agilent Technologies, Inc. 2009
Printed in USA, November 4, 2009
5989-6004EN RevA



Agilent Technologies